

# pump *in* style®

BREASTPUMP

*With Exclusive Medela 2-Phase  
Expression Technology®*

- Stylish microfiber shoulder bag
- Battery Pack (8AA batteries not included)
- Insulated PVC-free cooler carrier keeps milk cool for up to 12 hrs
- 4 collection containers & lids
- 2-24mm PersonalFit™ Breastshields
- 2 valves & 4 valve membranes
- Breastfeeding information & instruction guides



**CARE**  
MEDICAL &  
REHABILITATION  
EQUIPMENT



"Partners In Care"

**SEATTLE, WASHINGTON**  
4135 Stone Way North - 98103  
(206) 547-2200

**EVERETT, WASHINGTON**  
3010 Grand Avenue - 98201  
(425) 252-6456

**TACOMA, WASHINGTON**  
1902 Tacoma Avenue South - 98402  
(253) 474-1774

**OLYMPIA, WASHINGTON**  
3430 Pacific Avenue SE #A5 - 98501  
(360) 459-0999

**LONGVIEW, WASHINGTON**  
925 Vandercook Way - 98632  
(360) 636-3120

**VANCOUVER, WASHINGTON**  
1720 NE Andresen Road - 98661  
(360) 693-4777

**PORTLAND, OREGON**  
1877 NE 7th Avenue - 97212  
(503) 288-8174

**SALEM, OREGON**  
1944 McGilchrist SE - 97302  
(503) 378-1756

**EUGENE, OREGON**  
1280 Charnelton - 97401  
(541) 343-0304

**CARE**

MEDICAL &  
REHABILITATION  
EQUIPMENT

# Breastpumps

 **medela®**

#1 CHOICE OF HOSPITALS AND MOTHERS



[www.caremedical.com](http://www.caremedical.com)


*With the electric Pump In Style® breastpump in a stylish shoulder bag the active and working mother is fully equipped for pumping - any time, any place, every day. Ideal for frequent and extended use!*

## FEATURES & BENEFITS

- 2-Phase Expression Technology - pump automatically switches from Stimulation to Expression phase
- Feels more like your baby than any other retail pump
- One-touch “let-down” button designed for faster milk flow
- Customized pumping with one-knob adjustable speed and vacuum control
- Stylish bag contains everything you need to pump
- Battery Pack allows you to pump anywhere, anytime
- Used for single or double pumping
- All parts that come in contact with breastmilk are **BPA/DEHP Free!**

**pump *in* style®**  
BREASTPUMP

*2-Phase Expression Technology®  
Nature and Technology Perfectly Combined!*



### Stimulation Phase

Fast pumping rhythm starts milk flowing.  
(120 cycles per minute)

### Expression Phase

Slower pumping rhythm expresses milk gently  
& efficiently. (45-78 cycles per minute)

## 2-PHASE EXPRESSION

The new standard for unsurpassed naturalness, comfort & efficiency in electric breastpumps.

Breastmilk is best for the baby and the baby's natural sucking rhythm is best for the mother. In intensive pioneer work, Medela has converted the baby's intuitive knowledge into technological know-how.

Babies change the way they suckle at the mother's breast during breastfeeding by nursing in two phases: initially, they suckle fast and light to stimulate the milk ejection reflex and to start the milk flowing. Then, they switch to a slower and deeper suck to get as much milk as quickly as possible.

2-Phase Expression's advanced pumping pattern mimics this natural nursing pattern allowing for the expression of milk more gently and efficiently in less time.

